

Quality Control Log: High/Low Controls

Purpose: to verify and document performance of glucose meter, document maintenance

Month _____ Year _____

Unit _____ Meter serial: _____

Chemstrip bG lot #: _____ Strip code: _____ Expiration date: _____

Control Solutions – Acceptable Ranges

Low Control Solution
 Lot # _____
 Exp. Date _____
 Acceptable range:
 _____ to _____

High control solution
 Lot # _____
 Exp. Date _____
 Acceptable range:
 _____ to _____

Date	Low	High	Maintenance	Signature	Comment on any test out of control range
1					
2					
3					
4					
5					
6					
7					
8					
9					
10					
11					
12					
13					
14					
15					
16					
17					
18					
19					
20					
21					
22					
23					
24					
25					
26					
27					
28					
29					
30					
31					